

VAMP4 抗原（重组蛋白）

中文名称：VAMP4 抗原（重组蛋白）

英文名称：VAMP4 Antigen (Recombinant Protein)

别名：VAMP-4; VAMP24

储存：冷冻（-20℃）

相关类别：抗原

概述：

Full length fusion protein

技术规格：

Background:	Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. This protein may play a role in trans-Golgi network-to-endosome transport.
Tags:	His tag C-Terminus, GST tag N-Terminus
Expression system:	Escherichia coli
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Gene Accession:	BC007019
Swissprot:	O75379
Synonyms:	VAMP-4; VAMP24
Full name:	vesicle-associated membrane protein 4