

## WDR33 抗原（重组蛋白）

中文名称：WDR33 抗原（重组蛋白）

英文名称：WDR33 Antigen (Recombinant Protein)

别 名：NET14; WDC146

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human WDR33

技术规格：

<b>Full name:</b>	WD repeat domain 33
<b>Synonyms:</b>	NET14; WDC146
<b>Swissprot:</b>	Q9C0J8
<b>Gene Accession:</b>	BC013990
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartic acid (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene

regulation. This gene is highly expressed in testis and the protein is localized to the nucleus. This gene may play important roles in the mechanisms of cytodifferentiation and/or DNA recombination. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.