

兔抗 GPR157 多克隆抗体

中文名称: 兔抗 GPR157 多克隆抗体

英文名称: Anti-GPR157 rabbit polyclonal antibody

别名: G protein-coupled receptor 157

储存: 冷冻 (-20℃)

抗原: GPR157

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Orphan receptor that promotes neuronal differentiation of radial glial progenitors (RGPs). The activity of this receptor is mediated by a G(q)-protein that activates a phosphatidylinositol-calcium second messenger.
Applications:	ELISA, IHC
Name of antibody:	GPR157
Immunogen:	Synthetic peptide of human GPR157
Full name:	G protein-coupled receptor 157
SwissProt:	Q5UAW9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:

40-200



