

兔抗 GPR160 多克隆抗体

中文名称: 兔抗 GPR160 多克隆抗体

英文名称: Anti-GPR160 rabbit polyclonal antibody

别名: G protein-coupled receptor 160; GPCR1; GPCR150

储存: 冷冻 (-20℃)

抗原: GPR160

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Orphan receptor.
Applications:	ELISA, IHC
Name of antibody:	GPR160
Immunogen:	Synthetic peptide of human GPR160
Full name:	G protein-coupled receptor 160
Synonyms:	GPCR1; GPCR150
SwissProt:	Q9UJ42
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

