

ZNF471 抗原（重组蛋白）

中文名称：ZNF471 抗原（重组蛋白）

英文名称：ZNF471 Antigen (Recombinant Protein)

别名：ERP1; Z1971

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 427-626 amino acids of human ZNF471

技术规格：

Full name:	zinc finger protein 471
Synonyms:	ERP1; Z1971
Swissprot:	Q9BX82
Gene Accession:	BC125221
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in transcriptional regulation.