

## 兔抗 GPR180 多克隆抗体

中文名称: 兔抗 GPR180 多克隆抗体

英文名称: Anti-GPR180 rabbit polyclonal antibody

别名: G protein-coupled receptor 180; ITR

储存: 冷冻 (-20℃)

抗原: GPR180

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is a member of the G protein-coupled receptor superfamily. This protein is produced predominantly in vascular smooth muscle cells and may play an important role in the regulation of vascular remodeling.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GPR180
<b>Immunogen:</b>	Synthetic peptide of human GPR180
<b>Full name:</b>	G protein-coupled receptor 180
<b>Synonyms:</b>	ITR
<b>SwissProt:</b>	Q86V85

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	20-100



