

## ZNF581 抗原（重组蛋白）

中文名称：ZNF581 抗原（重组蛋白）

英文名称：ZNF581 Antigen (Recombinant Protein)

别 名：HSPC189

储 存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-197 amino acids of human ZNF581

技术规格：

|                           |                                                                                                 |
|---------------------------|-------------------------------------------------------------------------------------------------|
| <b>Full name:</b>         | zinc finger protein 581                                                                         |
| <b>Synonyms:</b>          | HSPC189                                                                                         |
| <b>Swissprot:</b>         | Q9P0T4                                                                                          |
| <b>Gene Accession:</b>    | BC010097                                                                                        |
| <b>Purity:</b>            | >85%, as determined by Coomassie blue stained SDS-PAGE                                          |
| <b>Expression system:</b> | Escherichia coli                                                                                |
| <b>Tags:</b>              | His tag C-Terminus, GST tag N-Terminus                                                          |
| <b>Background:</b>        | ZNF581 contains 4 C2H2-type zinc fingers. ZNF581 may be involved in transcriptional regulation. |