

ZNRD2 抗原(重组蛋白)

- 中文名称: ZNRD2 抗原(重组蛋白)
- 英文名称: ZNRD2 Antigen (Recombinant Protein)
- 别名: p27; SSSCA1
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述:

Fusion protein corresponding to a region derived from 2-199 amino acids of human ZNRD2

技术规格:

Full name:	zinc ribbon domain containing 2
Synonyms:	p27; SSSCA1
Swissprot:	O60232
Gene Accession:	BC014791
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This antigen is recognized by a subset of anti-centromere an tibodies from patients with scleroderma and/or Sjogren's syn drome. Subcellular localization has not yet been established. [provided by RefSeq, Jul 2008]