

兔抗 GRHL1 多克隆抗体

中文名称：兔抗 GRHL1 多克隆抗体

英文名称：Anti-GRHL1 rabbit polyclonal antibody

别名：grainyhead like transcription factor 1; MGR; NH32; LBP32; TFCP2L2

抗原：GRHL1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the grainyhead family of transcription factors. The encoded protein can exist as a homodimer or can form heterodimers with sister-of-mammalian grainyhead or brother-of-mammalian grainyhead. This protein functions as a transcription factor during development. [provided by RefSeq, Jun 2009]
Applications:	ELISA, WB, IHC
Name of antibody:	GRHL1
Immunogen:	Synthetic peptide of human GRHL1
Full name:	grainyhead like transcription factor 1

Synonyms:	MGR; NH32; LBP32; TFCP2L2
SwissProt:	Q9NZI5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	70 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	500-2000



