

## 兔抗 MANSC1 多克隆抗体

中文名称：兔抗 MANSC1 多克隆抗体

英文名称：Anti-MANSC1 rabbit polyclonal antibody

别名：MANSC domain containing 1; LOH12CR3; 9130403P13Rik

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MANSC1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	MANSC1 (MANSC domain-containing protein 1), also known as LOH12CR3 (Loss of heterozygosity 12 chromosomal region 3 protein), is a 414 amino acid single-pass membrane protein. Expressed throughout the body, MANSC1 contains one MANSC domain and is encoded by a gene that is located on chromosome 12. Encoding over 1,100 genes within 132 million bases, chromosome 12 makes up about 4.5% of the human genome. A number of skeletal deformities are linked to chromosome 12, including hypochondrogenesis, achondrogenesis and Kniest dysplasia. Chromosome 12 is also home to a ho
--------------------	---

	meobox gene cluster which encodes crucial transcription factors for morphogenesis, and the natural killer complex gene cluster encoding C-type lectin proteins which mediate the NK cell response to MHC I interaction.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MANSC1
<b>Immunogen:</b>	Fusion protein of human MANSC1
<b>Full name:</b>	MANSC domain containing 1
<b>Synonyms:</b>	LOH12CR3; 9130403P13Rik
<b>SwissProt:</b>	Q9H8J5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human cervical cancer
<b>IHC Recommend dilution:</b>	50-300



