

兔抗 MAP2 多克隆抗体

中文名称: 兔抗 MAP2 多克隆抗体

英文名称: Anti-MAP2 rabbit polyclonal antibody

别名: MAP2A; MAP2B; MAP2C

储存: 冷冻 (-20℃)

抗原: MAP2

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	MAP2

Immunogen:	Fusion protein of human MAP2
Full name:	microtubule associated protein 2
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
SwissProt:	P11137
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	200 kDa
WB Positive control:	Mouse brain tissue lysate
WB Recommended dilution:	1000-5000



