

兔抗 GRIPAP1 多克隆抗体

中文名称: 兔抗 GRIPAP1 多克隆抗体

英文名称: Anti-GRIPAP1 rabbit polyclonal antibody

别名: GRIP1 associated protein 1; GRASP-1

抗原: GRIPAP1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a guanine nucleotide exchange factor for the Ras family of small G proteins (RasGEF). The encoded protein interacts in a complex with glutamate receptor interacting protein 1 (GRIP1) and plays a role in the regulation of AMPA receptor function.
Applications:	ELISA, IHC
Name of antibody:	GRIPAP1
Immunogen:	Synthetic peptide of human GRIPAP1
Full name:	GRIP1 associated protein 1
Synonyms:	GRASP-1
SwissProt:	Q4V328

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300



