

## 兔抗 GRK1 多克隆抗体

中文名称：兔抗 GRK1 多克隆抗体

英文名称：Anti-GRK1 rabbit polyclonal antibody

别名：RK; RHOK; GPRK1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GRK1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates rhodopsin and initiates its deactivation. Defects in GRK1 are known to cause Oguchi disease 2 (also known as stationary night blindness Oguchi type-2).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GRK1
<b>Immunogen:</b>	Synthetic peptide of human GRK1
<b>Full name:</b>	G protein-coupled receptor kinase 1

<b>Synonyms :</b>	RK; RHOK; GPRK1
<b>SwissProt:</b>	Q15835
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200

