

## 兔抗 GRK6 多克隆抗体

中文名称：兔抗 GRK6 多克隆抗体

英文名称：Anti-GRK6 rabbit polyclonal antibody

别名：GPRK6

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：GRK6

宿主：Rabbit

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates the activated forms of G protein-coupled receptors thus initiating their deactivation. Several transcript variants encoding different isoforms have been described for this gene.
<b>Applications:</b>	WB, IF
<b>Name of antibody:</b>	GRK6
<b>Immunogen:</b>	Synthesized peptide derived from N-terminal of human GRK6.

<b>Full name:</b>	G protein-coupled receptor kinase 6
<b>Synonyms :</b>	GPRK6
<b>SwissProt:</b>	P43250
<b>WB Predicted band size:</b>	66 kDa
<b>WB Positive control:</b>	Jurkat cells, COLO205 cells and Huvec cells lysates
<b>WB Recommended dilution:</b>	500-3000
<b>IF positive control:</b>	Hela cells
<b>IF Recommend dilution:</b>	100-500



