

兔抗 GRWD1 多克隆抗体

中文名称：兔抗 GRWD1 多克隆抗体

英文名称：Anti-GRWD1 rabbit polyclonal antibody

别名：glutamate rich WD repeat containing 1; CDW4; GRWD; RRB1; WDR28

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GRWD1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a glutamate-rich protein that contains five WD-repeat motifs. The encoded protein may play a critical role in ribosome biogenesis and may also play a role in histone methylation through interactions with CUL4-DDB1 ubiquitin E3 ligase.
Applications:	ELISA, IHC
Name of antibody:	GRWD1
Immunogen:	Fusion protein of human GRWD1
Full name:	glutamate rich WD repeat containing 1
Synonyms:	CDW4; GRWD; RRB1; WDR28
SwissProt:	Q9BQ67

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-200

