

兔抗 GPR55 多克隆抗体

中文名称：兔抗 GPR55 多克隆抗体

英文名称：Anti-GPR55 rabbit polyclonal antibody

别名：G protein-coupled receptor 55; LPIR1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GPR55

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the G-protein-coupled receptor superfamily. The encoded integral membrane protein is a likely cannabinoid receptor. It may be involved in several physiological and pathological processes by activating a variety of signal transduction pathways.
Applications:	ELISA, IHC
Name of antibody:	GPR55
Immunogen:	Synthetic peptide of human GPR55
Full name:	G protein-coupled receptor 55
Synonyms:	LPIR1
SwissProt:	Q9Y2T6

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200





