

兔抗 GPR62 多克隆抗体

中文名称：兔抗 GPR62 多克隆抗体

英文名称：Anti-GPR62 rabbit polyclonal antibody

别名：G protein-coupled receptor 62; GPCR8; KPG_005

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GPR62

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Orphan receptor.Expressed in brain; detected in the basal forebrain, frontal cortex, caudate, putamen, thalamus and hippocampus.
Applications:	ELISA, WB, IHC
Name of antibody:	GPR62
Immunogen:	Synthetic peptide of human GPR62
Full name:	G protein-coupled receptor 62
Synonyms:	GPCR8; KPG_005
SwissProt:	Q9BZJ7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:	40-200
WB Predicted band size:	38 kDa
WB Positive control:	A549 cell lysate
WB Recommended dilution:	500-2000



