

兔抗 HCN3 多克隆抗体

中文名称：兔抗 HCN3 多克隆抗体

英文名称：Anti-HCN3 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：HCN3

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a multi-pass membrane protein that functions as a voltage gated cation channel. The encoded protein is a member of a family of closely related cyclic adenosine monophosphate-binding channel proteins. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	HCN3
Immunogen:	Synthetic peptide of human HCN3
Full name:	hyperpolarization activated cyclic nucleotide gated potassium channel 3
SwissProt:	Q9P1Z3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human tonsil

IHC Recommend dilution:	25-100
-------------------------	--------



