

兔抗 MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185)多克隆抗体

中文名称: 兔抗 MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185)多克隆抗体

英文名称: Anti-MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185) rabbit polyclonal

antibody

别 名: Stress-activated protein kinase JNK1; c-Jun N-terminal kinase 1; JNK-46

相关类别: 一抗

储 存: 冷冻 (-20℃) 避光

宿 主: Rabbit

抗 原: MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185)

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	ription factors, primarily components of AP-1 such as JUN,
	JDP2 and ATF2 and thus regulates AP-1 transcriptional activ
	ity. In T-cells, JNK1 and JNK2 are required for polarized diff
	erentiation of T-helper cells into Th1 cells By similarity. Pho
	sphorylates heat shock factor protein 4 (HSF4). /Responds t

o activation by environmental stress and pro-inflammatory

Responds to activation by environmental stress and pro-infl ammatory cytokines by phosphorylating a number of transc



	cytokines by phosphorylating a number of transcription fact ors, primarily components of AP-1 such as c-Jun and ATF2 and thus regulates AP-1 transcriptional activity. In T-cells, J NK1 and JNK2 are required for polarized differentiation of T-helper cells into Th1 cells. JNK2 isoforms display different binding patterns: alpha-1 and alpha-2 preferentially bind to c-Jun, whereas beta-1 and beta-2 bind to ATF2. However, t here is no correlation between binding and phosphorylation , which is achieved at about the same efficiency by all isof orms. JUNB is not a substrate for JNK2 alpha-2, and JUND binds only weakly to it.
Applications:	WB, IF
Name of antibody:	MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185)
Immunogen:	Synthetic peptide of human MAPK8/MAPK9/MAPK10 (phosp ho-Thr183/Tyr185)
Full name:	MAPK8/MAPK9/MAPK10 (phospho-Thr183/Tyr185)
Synonyms :	Stress-activated protein kinase JNK1 ; c-Jun N-terminal kina se 1; JNK-46
SwissProt:	P45983/P45984/P53779
WB Predicted band size:	46 kDa; 54 kDa
WB Positive control:	C6 cells treated with anisomycin
WB Recommended dilution:	500-1000
IF Positive control:	Hela cells
IF Recommended dilution	100-200





