

BTG3 抗原(重组蛋白)

中文名称: BTG3 抗原(重组蛋白)

英文名称: BTG3 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 22-252 amino acids of human BTG3

技术规格

Full name:	BTG anti-proliferation factor 3
Synonyms:	ANA; TOB5; TOFA; APRO4; TOB55; ANA/BTG3
Swissprot:	Q14201
Gene Accession:	BC011957
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear t o have antiproliferative properties. This encoded protein might play a role in neurogenesis in the central nervous system. Two t ranscript variants encoding different isoforms have been found f or this gene.