

兔抗 FRG1 多克隆抗体

中文名称：兔抗 FRG1 多克隆抗体

英文名称：Anti-FRG1 rabbit polyclonal antibody

别名：FSDH region gene 1; FSG1; FRG1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FRG1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene maps to a location 100 kb centromeric of the repeat units on chromosome 4q35 which are deleted in facioscapulohumeral muscular dystrophy (FSDH). It is evolutionarily conserved and has related sequences on multiple human chromosomes but DNA sequence analysis did not reveal any homology to known genes. In vivo studies demonstrate the encoded protein is localized to the nucleolus. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC
Name of antibody:	FRG1
Immunogen:	Fusion protein of human FRG1

Full name:	FSHD region gene 1
Synonyms:	FSG1; FRG1A
SwissProt:	Q14331
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human esophagus cancer
IHC Recommend dilution:	100-300



