

兔抗 MARVELD1 多克隆抗体

中文名称：兔抗 MARVELD1 多克隆抗体

英文名称：Anti-MARVELD1 rabbit polyclonal antibody

别名：GB14; MARVD1; MRVLDC1; bA548K23.8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MARVELD1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

MARVELD1 (MARVEL domain containing 1), also known as GB14, MARVD1 or MRVLDC1, is a 173 amino acid protein that is widely expressed and contains one MARVEL (MAL and related proteins for vesicle trafficking and membrane link) domain. The MARVEL domain, which is common in Occludin proteins, is a 130 amino acid motif that contains four transmembrane helices, both of which have cytoplasmic N- and C-terminal regions. Localizing to the nucleus, MARVELD1 is downregulated in primary multiple tumors of liver, bladder, ovary, vulva, breast, testis, uterus, kidney and testis. The gene encoding MARV

	ELD1 maps to human chromosome 10, which makes up approximately 4.5% of total DNA in cells and encodes nearly 1,200 genes. Defects in some of the genes that map to chromosome 10 are associated with Charcot-Marie Tooth disease, Jackson-Weiss syndrome, Usher syndrome, nonsyndromic deafness, Wolman's syndrome and Cowden syndrome.
Applications:	ELISA, WB, IHC
Name of antibody:	MARVELD1
Immunogen:	Synthetic peptide of human MARVELD1
Full name:	MARVEL domain containing 1
Synonyms :	GB14; MARVD1; MRVLDC1; bA548K23.8
SwissProt:	Q9BSK0
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human cervical cancer and Human colon cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	19 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	200-1000

