

兔抗 HES7 多克隆抗体

中文名称:	兔抗 HES7 多克隆抗体
英文名称:	Anti-HES7 rabbit polyclonal antibody
别名:	hes family bHLH transcription factor 7; SCDO4; bHLHb37
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	HES7
反应种属:	Human, Mouse
标记物:	Unconjugate
相关类别:	一抗
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes a member of the hairy and enhancer of split family of bHLH transcription factors. The mouse ortholog of this gene is regulated by Notch signaling. The protein functions as a transcriptional repressor, and is implicated in correct patterning of the axial skeleton. A mutation in this gene has been shown to result in spondylocostal dysostosis. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	HES7
Immunogen:	Synthetic peptide of human HES7

Full name:	hes family bHLH transcription factor 7
Synonyms:	SCDO4; bHLHb37
SwissProt:	Q9BYE0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human liver cancer
IHC Recommend dilution:	30-150



