

兔抗 MAST3 多克隆抗体

中文名称：兔抗 MAST3 多克隆抗体

英文名称：Anti-MAST3 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

抗 原：MAST3

宿 主：Rabbit

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family.
Applications:	WB
Name of antibody:	MAST3
Immunogen:	Synthesized peptide derived from internal of human MAST3 .
Full name:	microtubule associated serine/threonine kinase 3
Synonyms :	0
SwissProt:	O60307
WB Predicted band size:	143 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

