

CAP2 抗原（重组蛋白）

中文名称：CAP2 抗原（重组蛋白）

英文名称：CAP2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a C terminal 300 amino acids of human CAP2

技术规格

Full name:	cyclase associated actin cytoskeleton regulatory protein 2
Swissprot:	P40123
Gene Accession:	BC008481
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene was identified by its similarity to the gene for human adenyl cyclase-associated protein. The function of the protein encoded by this gene is unknown. However, the protein appears to be able to interact with adenyl cyclase-associated protein and actin.