

兔抗 HIPK2 多克隆抗体

中文名称: 兔抗 HIPK2 多克隆抗体

英文名称: Anti-HIPK2 rabbit polyclonal antibody

别名: PRO0593

抗原: HIPK2

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a conserved serine/threonine kinase that is a member of the homeodomain-interacting protein kinase family. The encoded protein interacts with homeodomain transcription factors and many other transcription factors such as p53, and can function as both a corepressor and a coactivator depending on the transcription factor and its subcellular localization. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	HIPK2
Immunogen:	Synthetic peptide of human HIPK2

Full name:	homeodomain interacting protein kinase 2
Synonyms:	PRO0593
SwissProt:	Q9H2X6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	25-100



