

兔抗 MCPH1 多克隆抗体

中文名称：兔抗 MCPH1 多克隆抗体

英文名称：Anti-MCPH1 rabbit polyclonal antibody

别名：MCT; BRIT1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：MCPH1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a DNA damage response protein. The encoded protein may play a role in G2/M checkpoint arrest via maintenance of inhibitory phosphorylation of cyclin-dependent kinase 1. Mutations in this gene have been associated with primary autosomal recessive microcephaly 1 and premature chromosome condensation syndrome. Alternatively spliced transcript variants have been described.
Applications:	WB, IHC
Name of antibody:	MCPH1
Immunogen:	Synthesized peptide derived from internal of human MCPH1.

Full name:	microcephalin 1
Synonyms :	MCT; BRIT1
SwissProt:	Q8NEM0
IHC positive control:	Human thymus glan tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	93 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000



