

DKKL1 抗原（重组蛋白）

中文名称：DKKL1 抗原（重组蛋白）

英文名称：DKKL1 Antigen (Recombinant Protein)

储 存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 31-242 amino acids of human DKKL1

技术规格

Full name:	dickkopf like acrosomal protein 1
Synonyms:	SGY; CT34; SGY1; SGY-1
Swissprot:	Q9UK85
Gene Accession:	BC030581
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The dickkopf protein family interacts with the Wnt signaling pathway and its members are characterized by two conserved cysteine-rich domains. This gene encodes a secreted protein that has low sequence similarity to the dickkopf-3 protein. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.