

兔抗 MGST3 多克隆抗体

中文名称：兔抗 MGST3 多克隆抗体

英文名称：Anti-MGST3 rabbit polyclonal antibody

别名：microsomal glutathione S-transferase 3; GST-III

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MGST3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the MAPEG (Membrane Associated Proteins in Eicosanoid and Glutathione metabolism) protein family. Members of this family are involved in the production of leukotrienes and prostaglandin E, important mediators of inflammation. This gene encodes an enzyme which catalyzes the conjugation of leukotriene A4 and reduced glutathione to produce leukotriene C4. This enzyme also demonstrates glutathione-dependent peroxidase activity towards lipid hydroperoxides.
Applications:	ELISA, IHC
Name of antibody:	MGST3

Immunogen:	Fusion protein of human MGST3
Full name:	microsomal glutathione S-transferase 3
Synonyms:	GST-III
SwissProt:	O14880
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	30-150

