

兔抗 MLH3 多克隆抗体

中文名称：兔抗 MLH3 多克隆抗体

英文名称：Anti-MLH3 rabbit polyclonal antibody

别 名：HNPCC7

抗 原：MLH3

储 存：冷冻（-20℃）避光

宿 主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:

This gene is a member of the MutL-homolog (MLH) family of DNA mismatch repair (MMR) genes. MLH genes are implicated in maintaining genomic integrity during DNA replication and after meiotic recombination. The protein encoded by this gene functions as a heterodimer with other family members. Somatic mutations in this gene frequently occur in tumors exhibiting microsatellite instability, and germline mutations have been linked to hereditary nonpolyposis colorectal cancer type 7 (HNPCC7). Several alternatively spliced transcript variants have been identified, but the full-length nature of only two transcript variants has been deter

	mined.
Applications:	WB
Name of antibody:	MLH3
Immunogen:	Synthesized peptide derived from internal of human MLH3.
Full name:	mutL homolog 3
Synonyms :	HNPCC7
SwissProt:	Q9UHC1
WB Predicted band size:	164 kDa
WB Positive control:	K562 cells and HT29 cells lysates
WB Recommended dilution:	500-3000

