

兔抗 MLXIPL 多克隆抗体

中文名称：兔抗 MLXIPL 多克隆抗体

英文名称：Anti-MLXIPL rabbit polyclonal antibody

级 别：MIO; CHREBP; MONDOB; WBSCR14; WS-bHLH; bHLHd14

储 存：冷冻（-20℃）

抗 原：MLXIPL

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a basic helix-loop-helix leucine zipper transcription factor of the Myc/Max/Mad superfamily. This protein forms a heterodimeric complex and binds and activates, in a glucose-dependent manner, carbohydrate response element (ChoRE) motifs in the promoters of triglyceride synthesis genes. The gene is deleted in Williams-Beuren syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at chromosome 7q11.23. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	MLXIPL
Immunogen:	Synthetic peptide of human MLXIPL
Full name:	MLX interacting protein-like
Synonyms:	MIO; CHREBP; MONDOB; WBSCR14; WS-bHLH; bHLHd14
SwissProt:	Q9NP71
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	20-100

