

兔抗 HMP19 多克隆抗体

中文名称: 兔抗 HMP19 多克隆抗体

英文名称: Anti-HMP19 rabbit polyclonal antibody

储 存: 冷冻 (-20℃) 避光

抗 原: HMP19

宿 主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Belongs to the NSG family.
Applications:	IHC
Name of antibody:	HMP19
Immunogen:	Synthesized peptide derived from internal of human NSG2.
Full name:	HMP19 protein
Synonyms :	0
SwissProt:	Q9Y328
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

