

## 兔抗 MMP15 多克隆抗体

中文名称：兔抗 MMP15 多克隆抗体

英文名称：Anti-MMP15 rabbit polyclonal antibody

别 名：MTMMP2, SMCP-2, MT2-MMP

抗 原：MMP15

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the protein encoded by this gene is a member of the membrane-type MMP (MT-MMP) subfamily; each member of this subfamily contains a potential transmembrane domain suggesting that these proteins are expressed at the
--------------------	---

	he cell surface rather than secreted.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	MMP15
<b>Immunogen:</b>	Fusion protein of human MMP15
<b>Full name:</b>	matrix metalloproteinase 15 (membrane-inserted)
<b>Synonyms :</b>	MTMMP2, SMCP-2, MT2-MMP
<b>SwissProt:</b>	P51511
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	76 kDa
<b>WB Positive control:</b>	Mouse stomach, heart and liver tissue, HT-29 cells
<b>WB Recommended dilution:</b>	500-2000

