

兔抗 GABRB1 多克隆抗体

中文名称：兔抗 GABRB1 多克隆抗体

英文名称：Anti-GABRB1 rabbit polyclonal antibody

别名：gamma-aminobutyric acid type A receptor beta1 subunit; EIEE45

相关类别：一抗

储存：冷冻（-20℃）

抗原：GABRB1

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The gamma-aminobutyric acid (GABA) A receptor is a multisubunit chloride channel that mediates the fastest inhibitory synaptic transmission in the central nervous system. This gene encodes GABA A receptor, beta 1 subunit. It is mapped to chromosome 4p12 in a cluster comprised of genes encoding alpha 4, alpha 2 and gamma 1 subunits of the GABA A receptor. Alteration of this gene is implicated in the pathogenetics of schizophrenia.
Applications:	ELISA, IHC
Name of antibody:	GABRB1
Immunogen:	Synthetic peptide of human GABRB1

Full name:	gamma-aminobutyric acid type A receptor beta1 subunit
Synonyms:	EIEE45
SwissProt:	P18505
ELISA Recommended dilution :	5000-10000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	50-200



