

兔抗 MMP19 多克隆抗体

中文名称：兔抗 MMP19 多克隆抗体

英文名称：Anti-MMP19 rabbit polyclonal antibody

别名：MMP18, RASI-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MMP19

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of proteins that are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded protein is secreted as an inactive proprotein, which is activated upon cleavage by extracellular proteases. Alternative splicing results in multiple transcript variants for this gene.
Applications:	ELISA, IHC
Name of antibody:	MMP19
Immunogen:	Synthetic peptide of human MMP19

Full name:	matrix metallopeptidase 19
Synonyms:	CODA; MMP18; RASI-1
SwissProt:	Q99542
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

