

## 兔抗 FOXS1 多克隆抗体

中文名称：兔抗 FOXS1 多克隆抗体

英文名称：Anti-FOXS1 rabbit polyclonal antibody

别名：forkhead box S1; FKHL18; FREAC10

抗原：FOXS1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FOXS1
<b>Immunogen:</b>	Fusion protein of human FOXS1
<b>Full name:</b>	forkhead box S1
<b>Synonyms:</b>	FKHL18; FREAC10
<b>SwissProt:</b>	O43638
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200



