

兔抗 MMP21 多克隆抗体

中文名称：兔抗 MMP21 多克隆抗体

英文名称：Anti-MMP21 rabbit polyclonal antibody

别名：no me

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MMP21

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Synthetic peptide of human MMP21
Applications:	ELISA, IHC
Name of antibody:	MMP21
Immunogen:	Synthetic peptide of human MMP21
Full name:	matrix metalloproteinase 21
Synonyms:	no me
SwissProt:	Q8N119
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer;Human thyroid cancer
IHC Recommend dilution:	25-100

