

## 兔抗 MMP24 多克隆抗体

中文名称：兔抗 MMP24 多克隆抗体

英文名称：Anti-MMP24 rabbit polyclonal antibody

别名：MMP25; MMP-24; MT5MMP; MTMMP5; MT-MMP5; MT5-MMP; MT-MMP 5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MMP24

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	MMP24
<b>Immunogen:</b>	Synthetic peptide of humanMMP24
<b>Full name:</b>	matrix metalloproteinase 24 (membrane-inserted)

<b>Synonyms :</b>	MMP25; MMP-24; MT5MMP; MTMMP5; MT-MMP5; MT5-MMP; MT-MMP 5
<b>SwissProt:</b>	Q9Y5R2
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	73 kDa
<b>WB Positive control:</b>	293T and HepG2 cells
<b>WB Recommended dilution:</b>	200-1000

