

兔抗 HNRNPA1L2 多克隆抗体

中文名称：兔抗 HNRNPA1L2 多克隆抗体

英文名称：Anti-HNRNPA1L2 rabbit polyclonal antibody

别名：heterogeneous nuclear ribonucleoprotein A1-like 2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HNRNPA1L2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Involved in the packaging of pre-mRNA into hnRNP particles, transport of poly(A) mRNA from the nucleus to the cytoplasm and may modulate splice site selection.
Applications:	ELISA, WB, IHC
Name of antibody:	HNRNPA1L2
Immunogen:	Synthetic peptide of human HNRNPA1L2
Full name:	heterogeneous nuclear ribonucleoprotein A1-like 2
SwissProt:	Q32P51
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human breast cancer

IHC Recommend dilution:	20-100
WB Predicted band size:	34 kDa
WB Positive control:	293T, K562, HepG2, 231, Hela and Jurkat cell lysates
WB Recommended dilution:	500-2000



