

兔抗 GSN 多克隆抗体

中文名称：兔抗 GSN 多克隆抗体

英文名称：Anti-GSN rabbit polyclonal antibody

别名：ADF, AGEL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GSN

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene binds to the "plus" ends of actin monomers and filaments to prevent monomer exchange. The encoded calcium-regulated protein functions in both assembly and disassembly of actin filaments. Defects in this gene are a cause of familial amyloidosis Finnish type (FAF). Multiple transcript variants encoding several different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	GSN
Immunogen:	Synthetic peptide of human GSN
Full name:	gelsolin

Synonyms :	ADF, AGEL
SwissProt:	P06396
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and Human thyroid cancer
IHC Recommend dilution:	25-100

