

兔抗 MMP26 多克隆抗体

中文名称：兔抗 MMP26 多克隆抗体

英文名称：Anti-MMP26 rabbit polyclonal antibody

别名：matrix metalloproteinase 26

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MMP26

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature enzyme. This enzyme may degrade collagen type IV, fibronectin, fibrinogen, and beta-casein, and activate matrix metalloproteinase-9 by cleavage. The protein differs from most MMP family members in that it lacks a conserved C-terminal protein domain. The enc

	oded protein may promote cell invasion in multiple human cancers.
Applications:	ELISA, IHC
Name of antibody:	MMP26
Immunogen:	Synthetic peptide of human MMP26
Full name:	matrix metalloproteinase 26
SwissProt:	Q9NRE1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human breast cancer
IHC Recommend dilution:	25-100



