

## 兔抗 MMP27 多克隆抗体

中文名称：兔抗 MMP27 多克隆抗体

英文名称：Anti-MMP27 rabbit polyclonal antibody

别 名：MMP-27

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：MMP27

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMPs are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MMP27
<b>Immunogen:</b>	Synthetic peptide of human MMP27
<b>Full name:</b>	Matrix metalloproteinase 27

<b>Synonyms :</b>	MMP-27
<b>SwissProt:</b>	Q9H306
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human breast cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

