

兔抗 HNRNPA3 多克隆抗体

中文名称: 兔抗 HNRNPA3 多克隆抗体

英文名称: Anti-HNRNPA3 rabbit polyclonal antibody

别名: heterogeneous nuclear ribonucleoprotein A3; FBRNP; HNRPA3; D10S102; 2610510D13Rik

抗原: HNRNPA3

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Plays a role in cytoplasmic trafficking of RNA. Binds to the c is-acting response element, A2RE. May be involved in pre-m RNA splicing.
Applications:	ELISA, IHC
Name of antibody:	HNRNPA3
Immunogen:	Synthetic peptide of human HNRNPA3
Full name:	heterogeneous nuclear ribonucleoprotein A3
Synonyms:	FBRNP; HNRPA3; D10S102; 2610510D13Rik
SwissProt:	P51991
ELISA Recommended diluti	5000-10000

on:	
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	25-100



