

兔抗 GSPT2 多克隆抗体

中文名称：兔抗 GSPT2 多克隆抗体

英文名称：Anti-GSPT2 rabbit polyclonal antibody

别名：G1 to S phase transition 2; GST2; ERF3B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GSPT2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a GTPase that belongs to the GTP-binding elongation factor family. The encoded protein is a polypeptide release factor that complexes with eukaryotic peptide chain release factor 1 to mediate translation termination. This protein may also be involved in mRNA stability.
Applications:	ELISA, WB, IHC
Name of antibody:	GSPT2
Immunogen:	Fusion protein of human GSPT2
Full name:	G1 to S phase transition 2
Synonyms:	GST2; ERF3B

SwissProt:	Q8IYD1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	40-200
WB Predicted band size:	69 kDa
WB Positive control:	PC3 cell lysate
WB Recommended dilution:	500-2000

kDa

250—

130—

95—

72—

55—

36—

28—



