

兔抗 HNRNPD(Ab-83) 多克隆抗体

- 中文名称：兔抗 HNRNPD(Ab-83) 多克隆抗体
- 英文名称：Anti-HNRNPD(Ab-83) rabbit polyclonal antibody
- 别名：P37; AUF1; AUF1A; HNRPD; hnRNP0
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：HNRNPD(Ab-83)
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are nucleic acid binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport.
Applications:	WB, IHC, IF
Name of antibody:	HNRNPD(Ab-83)
Immunogen:	Synthesized non-phosphopeptide derived from human hnRNP around the phosphorylation site of serine 83 (N-S-S(p)-P-R).

Full name:	heterogeneous nuclear ribonucleoprotein D (AU-rich element R NA binding protein 1, 37kDa)
Synonyms :	P37; AUF1; AUF1A; HNRPD; hnRNP D0
SwissProt:	Q14103
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution :	50-100
WB Predicted band size:	38 kDa
WB Positive control:	293, Jurkat, 3T3, COLO cell lysates
WB Recommended diluti on:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



