

兔抗 HNRNPR 多克隆抗体

中文名称：兔抗 HNRNPR 多克隆抗体

英文名称：Anti-HNRNPR rabbit polyclonal antibody

别名：heterogeneous nuclear ribonucleoprotein R; HNRPR; hnRNP-R

相关类别：一抗

储存：冷冻（-20℃）

抗原：HNRNPR

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10.
Applications:	ELISA, IHC
Name of antibody:	HNRNPR
Immunogen:	Synthetic peptide of human HNRNPR
Full name:	heterogeneous nuclear ribonucleoprotein R
Synonyms:	HNRPR; hnRNP-R
SwissProt:	O43390

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	25-100



