

兔抗 GALC 多克隆抗体

中文名称：兔抗 GALC 多克隆抗体

英文名称：Anti-GALC rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：GALC

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a lysosomal protein which hydrolyzes the galactose ester bonds of galactosylceramide, galactosylsphingosine, lactosylceramide, and monogalactosyldiglyceride. Mutations in this gene have been associated with Krabbe disease, also known as globoid cell leukodystrophy. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Applications:	ELISA, WB, IHC
Name of antibody:	GALC
Immunogen:	Synthetic peptide of human GALC
Full name:	galactosylceramidase
SwissProt:	P54803

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	77 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000



