

兔抗 GALP 多克隆抗体

中文名称：兔抗 GALP 多克隆抗体

英文名称：Anti-GALP rabbit polyclonal antibody

别名：galanin like peptide

相关类别：一抗

抗原：GALP

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the galanin family of neuropeptides. The encoded protein binds galanin receptors 1, 2 and 3 with the highest affinity for galanin receptor 3 and has been implicated in biological processes involving the central nervous system including hypothalamic regulation of metabolism and reproduction. A peptide encoded by a splice variant of this gene, termed alarin, has vasoactive properties, displays antimicrobial activity against E. coli, and may serve as a marker for neuroblastic tumors.
Applications:	ELISA, WB, IHC

Name of antibody:	GALP
Immunogen:	Synthetic peptide of human GALP
Full name:	galanin like peptide
SwissProt:	Q9UBC7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	13 kDa
WB Positive control:	Human fetal brain tissue and Mouse thymus tissue lysates
WB Recommended dilution:	500-2000



