

兔抗 GTPBP10 多克隆抗体

中文名称：兔抗 GTPBP10 多克隆抗体

英文名称： Anti-GTPBP10 rabbit polyclonal antibody

别名： ObgH2; UG0751c10

相关类别： 一抗

抗原： GTPBP10

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Small G proteins, such as GTPBP10, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction
Applications:	ELISA, WB
Name of antibody:	GTPBP10
Immunogen:	Fusion protein of human GTPBP10
Full name:	GTP-binding protein 10
Synonyms :	ObgH2; UG0751c10
SwissProt:	A4D1E9

ELISA Recommended dilution:	1000-2000
WB Predicted band size:	43 kDa
WB Positive control:	A549 cells
WB Recommended dilution:	200-1000

