

兔抗 GUCA2A 多克隆抗体

- 中文名称：兔抗 GUCA2A 多克隆抗体
- 英文名称：Anti-GUCA2A rabbit polyclonal antibody
- 别名：guanylate cyclase activator 2A; GUCA2; STARA; GCAP-I
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：GUCA2A
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Endogenous activator of intestinal guanylate cyclase. It stimulates this enzyme through the same receptor binding region as the heat-stable enterotoxins.
Applications:	ELISA, WB, IHC
Name of antibody:	GUCA2A
Immunogen:	Synthetic peptide of human GUCA2A
Full name:	guanylate cyclase activator 2A
Synonyms:	GUCA2; STARA; GCAP-I
SwissProt:	Q02747
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer

IHC Recommend dilution:	30-150
WB Predicted band size:	12 kDa
WB Positive control:	Mouse small intestines tissue and K562 cell lysates
WB Recommended dilution:	200-1000



